

ITPRIPL2 Antibody, Biotin conjugated



PACO64286

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ITPRIPL2

Uniprot

Q3MIP1

Synonyms:

Inositol 1,4,5-trisphosphate receptor-interacting protein-like 2, ITPRIPL2

Immunogen:

Recombinant Human Inositol 1,4,5-trisphosphate receptor-interacting protein-like 2 protein (75-288AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A