

# Mmp2 Antibody, FITC conjugated



PACO64305

---

## Product Information

**Size:**

50ul

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Mmp2

**Uniprot**

P33434

**Synonyms:**

72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) [Cleaved into: PEX], Mmp2

**Immunogen:**

Recombinant Mouse 72 kDa type IV collagenase protein (412-662AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A