

OR10P1 Antibody



PACO64363

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

Odorant receptor.

Gene ID:

OR10P1

Uniprot

Q8NGE3

Synonyms:

Olfactory receptor 10P1 (Olfactory receptor 10P2) (Olfactory receptor 10P3) (Olfactory receptor OR12-7), OR10P1, OR10P1P OR10P2P OR10P3P


Immunogen:

Peptide sequence from Human Olfactory receptor 10P1 protein (260-272AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: HL60 whole cell lysate. All lanes: OR10P1 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 35 kDa. Observed band size: 35 kDa.