

OR2S2 Antibody



PACO64383

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:100-1:500

Protein Background:

Odorant receptor.

Gene ID:

OR2S2

Uniprot

Q9NQN1

Synonyms:

Olfactory receptor 2S2 (Olfactory receptor OR9-3), OR2S2


Immunogen:

Peptide sequence from Human Olfactory receptor 2S2 protein (262-276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescence staining of HeLa cells with PACO64383 at 1:200, counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).