

OR4D10 Antibody



PACO64395

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

Odorant receptor.

Gene ID:

OR4D10

Uniprot

Q8NGI6

Synonyms:

Olfactory receptor 4D10 (Olfactory receptor OR11-251), OR4D10, OR4D10P


Immunogen:

Peptide sequence from Human Olfactory receptor 4D10 protein (86-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: HepG2 whole cell lysate, A549 whole cell lysate, NTERA-2 whole cell lysate. All lanes: OR4D10 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 36 kDa. Observed band size: 36 kDa.