

OR5R1 Antibody



PACO64411

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

Odorant receptor.

Gene ID:

OR5R1

Uniprot

Q8NH85

Synonyms:

Olfactory receptor 5R1 (Olfactory receptor OR11-185), OR5R1, OR5R1P


Immunogen:

Peptide sequence from Human Olfactory receptor 5R1 protein (292-307AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: K562 whole cell lysate. All lanes: OR5R1 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 37 kDa. Observed band size: 37 kDa.