

OR5T2 Antibody, HRP conjugated



PACO64416

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Odorant receptor.

Gene ID:

OR5T2

Uniprot

Q8NKG2

Synonyms:

Olfactory receptor 5T2 (Olfactory receptor OR11-177), OR5T2

Immunogen:

Peptide sequence from Human Olfactory receptor 5T2 protein (298-310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A