

OR5T2 Antibody, Biotin conjugated



PACO64418

Product Information

Size:	Protein Background:
50ul	Odorant receptor.
Reactivity:	Gene ID:
Human	OR5T2
Source:	Uniprot
Rabbit	Q8NKG2
Isotype:	Synonyms:
IgG	Olfactory receptor 5T2 (Olfactory receptor OR11-177), OR5T2
Applications:	Immunogen:
ELISA	Peptide sequence from Human Olfactory receptor 5T2 protein (298-310AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A