

PCDHA10 Antibody, FITC conjugated



PACO64441

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

Gene ID:

PCDHA10

Uniprot

Q9Y5I2

Synonyms:

Protocadherin alpha-10 (PCDH-alpha-10), PCDHA10, CNRS8

Immunogen:

Recombinant Human Protocadherin alpha-10 protein (276-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A