

# PCDHGB7 Antibody, Biotin conjugated



PACO64454

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

**Gene ID:**

PCDHGB7

**Uniprot**

Q9Y5F8

**Synonyms:**

Protocadherin gamma-B7 (PCDH-gamma-B7), PCDHGB7

**Immunogen:**

Recombinant Human Protocadherin gamma-B7 protein (719-929AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A