

SUMF2 Antibody, HRP conjugated



PACO64536

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Lacks formyl-glycine generating activity and is unable to convert newly synthesized inactive sulfatases to their active form. Inhibits the activation of sulfatases by SUMF1.

Gene ID:

SUMF2

Uniprot

Q8NBJ7

Synonyms:

Sulfatase-modifying factor 2 (C-alpha-formylglycine-generating enzyme 2), SUMF2

Immunogen:

Recombinant Human Sulfatase-modifying factor 2 protein (26-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A