

# TAS2R7 Antibody, HRP conjugated



PACO64544

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Gustducin-coupled receptor implicated in the perception of bitter compounds in the oral cavity and the gastrointestinal tract. Signals through PLCB2 and the calcium-regulated cation channel TRPM5.

**Gene ID:**

TAS2R7

**Uniprot**

Q9NYW3

**Synonyms:**

Taste receptor type 2 member 7 (T2R7) (Taste receptor family B member 4) (TRB4), TAS2R7

**Immunogen:**

Peptide sequence from Human Taste receptor type 2 member 7 protein (151-168AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A