

# TBC1D3G Antibody



PACO64547

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:100-1:500

**Protein Background:**

Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.

**Gene ID:**

TBC1D3G

**Uniprot**

Q6DHY5

**Synonyms:**

TBC1 domain family member 3G, TBC1D3G

**Immunogen:**


Recombinant Human TBC1 domain family member 3G protein (312-549AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunofluorescence staining of U251 cells with PACO64547 at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).