

TBC1D3G Antibody



PACO64551

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:100-1:500

Protein Background:

Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.

Gene ID:

TBC1D3G

Uniprot

Q6DHY5

Synonyms:

TBC1 domain family member 3G, TBC1D3G


Immunogen:

Recombinant Human TBC1 domain family member 3G protein (312-549AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescence staining of U251 cells with PACO64551 at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).