



# Recombinant Protein Technical Manual

## Recombinant Human NRN1L Protein (His Tag)

RPES0240

### Product Data:

**Product SKU:** RPES0240

**Size:** 10µg

**Species:** Human

**Expression host:** Human Cells

**Uniprot:** Q496H8

### Protein Information:

**Molecular Mass:** 12.3 kDa

**AP Molecular Mass:** 14&16 kDa

**Tag:** C-6His

**Bio-activity:**

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Endotoxin:** < 1.0 EU per µg as determined by the LAL method.

**Storage:** Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Formulation:** Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

**Reconstitution:** Please refer to the printed manual for detailed information.

**Application:**

**Synonyms:** Neuritin-like protein; NRN1L; UNQ2446/PRO5725

## Immunogen Information:

**Sequence:** Ala36-Ala139

## Background:

Neuritin-like protein belongs to the neuritin family. Neuritin is a GPI-anchored protein that promotes neurite outgrowth and branching of neuritic processes in primary hippocampal and cortical cells. Neuritin expression also enhances the development of motor neuron axon arbors by promoting neuromuscular synaptogenesis and by stimulating the addition of new axon branches. Neuritin is induced by neuronal activity and by the neurotrophins, BDNF and NT3. NRN1L contains a consensus cleavage signal found in glycosylphosphatidylinositol (GPI)-anchored proteins.