



Recombinant Protein Technical Manual

Recombinant Human PPY Protein (His Tag)

RPES1172

Product Data:

Product SKU: RPES1172

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: P01298

Protein Information:

Molecular Mass: 7.8 kDa

AP Molecular Mass: 124 kDa

Tag: C-6His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Pancreatic prohormone;Pancreatic polypeptide;PNP;PPY

Immunogen Information:

Sequence: Ala30-Arg88

Background:

PPY belongs to the NPY family and is synthesized as a 95 aa polypeptide precursor in the pancreatic islets of Langerhans. It is cleaved into two peptide products; the active hormone of 36 aa and an icosapeptide of unknown function. The hormone acts as a regulator of pancreatic and gastrointestinal functions and may be important in the regulation of food intake. Plasma level of this hormone has been shown to be reduced in conditions associated with increased food intake and elevated in anorexia nervosa.