



Recombinant Protein Technical Manual
Recombinant Mouse Serpinb6b Protein (His Tag)
RPES2619

Product Data:

Product SKU: RPES2619

Size: 20µg

Species: Mouse

Expression host: Baculovirus-Insect Cells

Uniprot: NP_035584.1

Protein Information:

Molecular Mass: 44 kDa

AP Molecular Mass: 43 kDa

Tag: C-His

Bio-activity:

Purity: > 92 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 50mM Tris, 100mM NaCl, pH 7.0, 10% gly

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: NK13;ovalbumin;Spi12

Immunogen Information:

Sequence: Met 1-Pro 377

Background: