



Recombinant Protein Technical Manual
Recombinant Human KIAA1279 Protein (His & GST
Tag)
RPES3354

Product Data:

Product SKU: RPES3354

Size: 20µg

Species: Human

Expression host: Baculovirus-Insect Cells

Uniprot: Q96EK5

Protein Information:

Molecular Mass: 99.6 kDa

AP Molecular Mass: 9202 kDa

Tag: N-His & GST

Bio-activity:

Purity: > 90 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: KBP;KIAA1279;TTC20

Immunogen Information:

Sequence: Met 1-Thr621

Background: