



Recombinant Protein Technical Manual
Recombinant Human BTNL9/BTN3 Protein (His Tag)
RPES4051

Product Data:

Product SKU: RPES4051

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: Q6UXG8

Protein Information:

Molecular Mass: 15 kDa

AP Molecular Mass: 18-30 kDa

Tag: C-His

Bio-activity:

Purity: > 95% as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Butyrophilin-like protein 9; BTNL9

Immunogen Information:

Sequence: Glu48-Ser160

Background:

Butyrophilin-Like Protein 9 (BTNL9) is single-pass type I membrane protein member of the BTN/MOG family that belongs to the immunoglobulin superfamily. BTNL9 consists of two domains: one B30.2/SPRY domain and one Ig-like V-type (immunoglobulin-like) domain. Human BTNL9 mRNA has been identified in adipose, lung, thymus, spleen, colon, and cardiac tissues, but its highest levels of expression were found in B cells. BTNL9 expression has also been found to be down-regulated in colon cancer tumors.