



Recombinant Protein Technical Manual

Recombinant Human TMEM25 Protein (His Tag)

RPES4220

Product Data:

Product SKU: RPES4220

Size: 20µg

Species: Human

Expression host: HEK293 Cells

Uniprot: AAH51841.1

Protein Information:

Molecular Mass: 22.8 kDa

AP Molecular Mass: 43 kDa

Tag: C-His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: bA83N9.1

Immunogen Information:

Sequence: Met 1-Pro224

Background:

TMEM25 is a novel member of the immunoglobulin superfamily. Immunoglobulin superfamily members are implicated in immune responses, growth factor signaling, and cell adhesion. TMEM25 contains 1 Ig-like (immunoglobulin-like) domain and is a target of pharmacogenomics in the field of oncology and regenerative medicine. TMEM25 isoform 1, consisting of exons 1-9, encoded a 366-aa transmembrane protein. TMEM25 isoform 2, consisting of exons 1-4 and 6-9, encoded a 322-aa secreted protein. TMEM25 mRNA was expressed in brain, including cerebellar cortex and hippocampus, as well as in neuroblastoma, brain tumors, and gastric cancer. Human TMEM25 gene was located at the 11q23.3 oncogenomic recombination hotspot around the MLL amplicon and the neuroblastoma deleted region.