



Recombinant Protein Technical Manual

Recombinant Human C12orf53 Protein (Fc Tag)

RPES4251

Product Data:

Product SKU: RPES4251

Size: 20µg

Species: Human

Expression host: HEK293 Cells

Uniprot: Q8IYJ0

Protein Information:

Molecular Mass: 42.3 kDa

AP Molecular Mass: 53 kDa

Tag: C-Fc

Bio-activity:

Purity: > 96 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: C12orf53;leda;LEDA1;PANP;UNQ828/PRO1755

Immunogen Information:

Sequence: Met 1-Pro178

Background:

C12orf53 is mainly expressed in adult brain and cerebellum. It also can be detected in fetal brain and virtually no expression in spleen, heart, kidney, liver and dorsal ganglion relative to brain. C12orf53 acts as a ligand for PILRA in neural tissues, where it may be involved in immune regulation. Chromosome 12 encodes over 1,100 genes within 132 million bases. A number of skeletal deformities are linked to chromosome 12 including hypochondrogenesis, achondrogenesis and Kniest dysplasia. Noonan syndrome, which includes heart and facial developmental defects among the primary symptoms, is caused by a mutant form of PTPN11 gene product, SH-PTP2. Chromosome 12 is also home to a homeobox gene cluster which encodes crucial transcription factors for morphogenesis, and the natural killer complex gene cluster encoding C-type lectin proteins which mediate the NK cell response to MHC class I interaction.