

# Zea Mays CK2a Recombinant Protein



RPPB2690

---

## Product Information    Protein Information

**Product SKU:**

RPPB2690

**Accession:**

P28523

**Host:**

Escherichia Coli.

**Protein description:**

Casein Kinase 2 alpha Zea Mays Recombinant is a non-glycosylated polypeptide having a molecular mass of 39.2 kDa. Casein Kinase 2 is purified by proprietary chromatographic techniques.

**Synonyms:**

Casein kinase II subunit alpha, EC 2.7.11.1, CK II, CK2-alpha, CK2?.

**Formulation:**

CK2a is supplied in 50% glycerol.

**Purity:**

Greater than 99% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.