

Human C4c Recombinant Protein



RPPB3095

Product Information Protein Information

Product SKU:

RPPB3095

Protein description:

Human C4c produced in Human plasma having a molecular mass of 210 kDa.

Host:

Human Plasma.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Complement C4c, Complement Component C4c, C4c.

Formulation:

The Human Complement C4c solution is in Sodium Phosphate buffer pH 7.0 containing 0.15M NaCl and 0.09% NaN₃.

Purity:

Greater than 99.0%.

Stability:

Human C4c although stable at 4°C for 1 week, should be stored at -15°C. Please avoid freeze-thaw cycles.