

Human Cadherin-E Recombinant Protein



RPPB3124

Product Information Protein Information

Product SKU:

RPPB3124

Accession:

P12830

Host:

Escherichia Coli.

Protein description:

Cadherin-E Human Recombinant (aa 600-707) expressed in E.coli, shows a 38 kDa band on SDS-PAGE. The Cadherin-E is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Epithelial cadherin, E-cadherin, Uvomorulin, Cadherin-1, CAM 120/80, CD324 antigen, CDH1, CDHE, UVO, ECAD, LCAM, Arc-1, CD324, Cadherin-E.

Formulation:

Cadherin-E protein (100µg/1ml) in 50mM Tris-HCl, pH7.5 and 10mM L-glutathione (reduced).

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.

Amino Acid Sequence:

iffcernpkpq viniidadlp pntspftael thgasanwti qyndptqesi ilkpkmalev gdykinlklm dnqkdkvtt levsdcceg aagvcrkaqp veaglqi.