

# Human KRT8 Recombinant Protein



RPPB3824

---

## Product Information    Protein Information

**Product SKU:**

RPPB3824

**Accession:**

P05787

**Host:**

Escherichia Coli.

**Protein description:**

Cytokeratin 8 Human Recombinant full length protein expressed in E.coli, shows a 78 kDa SDS-PAGE. The Cytokeratin 8 is fused to GST-Tag and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Keratin type II cytoskeletal 8, Cytokeratin-8, CK-8, Keratin-8, K8, KRT8, CYK8, KO, CK8, K2C8, CARD2.

**Formulation:**

Cytokeratin-8 in 50mM Tris-HCl pH-7.5, 10mM L-glutathione (reduced).

**Stability:**

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.