

Human HbA1c Recombinant Protein



RPPB3827

Product Information Protein Information

Product SKU:

RPPB3827

Protein description:

The Human Hemoglobin A1c was purified from Human Erythrocytes.

Host:

Human Erythrocytes.

Appearance:

Clear red frozen solution.

Formulation:

The HbA1c is supplied in a proprietary buffer formulation having pH-8.0. The HbA1c was dispensed on the basis of HbA1c not total hemoglobin.

Purity:

Greater than 96.0%.

Stability:

Human HbA1c although stable at 4°C for 1 week, should be stored at -20°C.