

# Human LA/SS-B Recombinant Protein



RPPB3830

---

## Product Information    Protein Information

**Product SKU:**

RPPB3830

**Accession:**

P05455

**Host:**

Sf9 insect cells.

**Protein description:**

LA/SS-B Human Recombinant produced in SF9 is a single, glycosylated, polypeptide chain having a calculated molecular mass of 48 kDa. The LA/SS-B is expressed with a -6x His tag at N-terminus and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Lupus La protein, Sjogren syndrome type B antigen, SS-B, La ribonucleoprotein, La autoantigen, SSB, La, LARP3, LA/SS-B, La(SS-B).

**Formulation:**

The protein solution contains 20mM HEPES, pH 7.5, 400mM NaCl, 20% Glycerol.

**Purity:**

Greater than 80.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.