

Human HDL Recombinant Protein



RPPB3835

Product Information Protein Information

Product SKU:

RPPB3835

Protein description:

Human High Density Lipoprotein (HDL) produced in Human plasma.

Host:

Human plasma.

Appearance:

Yellow to orange liquid.

Synonyms:

High Density Lipoprotein, HDL.

Purity:

Cholesterol level >200mg/l. <1% LDL.

Stability:

Human HDL although stable at 4°C for 1 week, should be stored below -15°C (short term i.e. < 3 months) and below -70°C for long term.