

Human P Selectin Recombinant Protein



RPPB4317

Product Information Protein Information

Product SKU:

RPPB4317

Accession:

P16109

Host:

Escherichia Coli.

Protein description:

P-Selectin Human Recombinant is expressed in E. coli containing 566 amino acids 197-761 fused to an amino terminal hexahistidine tag.

Appearance:

Sterile Filtered liquid formulation.

Synonyms:

P-selectin, Granule membrane protein 140, GMP-140, PADGEM, Leukocyte-endothelial cell adhesion molecule 3, LECAM3, CD62 antigen-like family member P, CD62P antigen, SELP, GMRP, GRMP, CD62, PSEL, CD62P, GMP140, FLJ45155.

Formulation:

P-Selectin is supplied in 50% Glycerol.

Purity:

Greater than 95% as determined by SDS-PAGE. Single band on Western Blot.

Stability:

Store at 4°C if entire vial will be used within 1-2 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.