

Human RPLP0 Recombinant Protein



RPPB4501

Product Information Protein Information

Product SKU:

RPPB4501

Accession:

P05388

Host:

Sf9 insect cells.

Protein description:

Ribosomal Phosphoprotein P0 Human Recombinant produced in SF9 is a glycosylated, polypeptide chain having a molecular mass of 35,096 Dalton. RPLP0 is expressed with a -6xHis tag and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

60S acidic ribosomal protein P0, Ribosomal Phosphoprotein P0, L10E, RPLP0, Ribosomal Protein Large P0, RPP0, P0, PRLP0, MGC88175, MGC111226.

Formulation:

RPLP0 is supplied in 20mM HEPES buffer pH-7.5, 0.01mM EDTA & 0.02% SDS.

Purity:

Greater than 95% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.