

Human RPLP1 Recombinant Protein



RPPB4502

Product Information Protein Information

Product SKU:

RPPB4502

Accession:

P05386

Host:

Sf9 insect cells.

Protein description:

RPLP1 is a full-length cDNA coding for the human ribosomal P1 phosphoprotein having a molecular mass of 12,336 Dalton (pH 4.75). RPLP1 protein is fused to a hexa-histidine purification tag.

Synonyms:

RPP1, RRP1, 60S acidic ribosomal protein P1, FLJ27448, MGC5215, acidic ribosomal phosphoprotein P1.

Formulation:

RPLP1 (0.34mg/ml) is supplied in 20mM HEPES buffer pH-8.0, 200mM NaCl and 6M Urea.

Purity:

Greater than 80% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.