

Human SNCA 61-140 Recombinant Protein



RPPB4730

Product Information Protein Information

Product SKU:

RPPB4730

Accession:

P37840

Host:

Escherichia Coli.

Protein description:

A-Synuclein 61-140 Human Recombinant which is a deletion mutant of the a-synuclein amino acids 61-140, produced in E.Coli is a single, non-glycosylated polypeptide chain of 81 amino acids having a molecular mass of 8.4kDa (molecular size on SDS-PAGE will appear higher), with an additional Met attached at the N-terminus. The Recombinant Human a-Synuclein 61-140 is purified by proprietary chromatographic techniques.

Appearance:

Sterile filtered colorless solution.

Synonyms:

Alpha-synuclein, Non-A beta component of AD amyloid, Non-A4 component of amyloid precursor, NACP, PD1, PARK1, PARK4, MGC110988, a-Synuclein, SNCA.

Formulation:

The SNCA 61-140 protein solution (1mg/ml) contains 20mM Tris-HCl buffer pH 7.5 and 100mM NaCl.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MEQVTNVGGAV VTGVTAVAQK TVEGAGSIAA ATGFVKKDKL GKNEEGAPQE GILEDMPVDP DNEAYEMPSE EGYQDYEPEA.