

TBEV NE Recombinant Protein



RPPB5626

Product Information Protein Information

Product SKU:

RPPB5626

Protein description:

The E.coli derived recombinant 37 kDa protein NE contains the Tick-borne encephalitis virus N-terminus regions of glycoprotein E.

Formulation:

20mM MES pH 6.5, 8M urea, 200 mM NaCl and 0.05% Tween-20.

Purity:

Encephalitis protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Encephalitis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.