

TriMethyl-Histone H3-K4 Monoclonal Antibody

CAB22146



Product Information

Product SKU:

CAB22146

Gene ID:

82908350

Observed MW:

17kDa

Calculated MW:

15kDa

Category:

Primary Antibody

Uniprot:

Q16695P68431

Reactivity:Human, Mouse, Rat, Other
(Wide Range Predicted)**Purification Method:**

Affinity purification

Source:

Rabbit

Isotype:

IgG

Immunogen information

Background:

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is located separately from the other H3 genes that are in the histone gene cluster on chromosome 6p22-p21.3.

Size:

20uL, 100uL

Applications:

Western blotting Immunohistochemistry Immunofluorescence

Synonyms:

H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; TriMethyl-Histone H3-K4

Immunogen:

A synthetic trimethylated peptide around K4 of human Histone H3 (NP_003520.1).

Recommended Dilution:DB 1:10000 - 1:60000 WB 1:1000- 1:5000 IHC-P 1:50 - 1:200 IF/ICC 1:50 - 1:200 ChIP 5µg antibody for 5µg-10µg of ChromatinCUT&Tag 10⁵ cells /1 µgChIP-seq 1:50 - 1:200**Storage:**

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Contact Details | Dublin, Ireland**Email:** techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.