

**HDFP984**

---

## Product Information

<b>Product SKU:</b>	HDFP984	<b>Expression Host:</b>	HEK293	<b>Size:</b>	10µg
<b>Target:</b>	GHRHR	<b>Tag:</b>	C-Flag&Strep Tag		

---

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot ID:</b>	Q02643
<b>Molecular Weight:</b>	The human full length GHRHR-Strep protein has a MW of 47.4 kDa		

---

## Protein Information

**Background:** This gene encodes a receptor for growth hormone-releasing hormone. Binding of this hormone to the receptor leads to synthesis and release of growth hormone. Mutations in this gene have been associated with isolated growth hormone deficiency (IGHD), also known as Dwarfism of Sindh, a disorder characterized by short stature. [provided by RefSeq, Jun 2010]

**Synonyms:** GHRFR, GRFR, IGHD1B, IGHD4

**Protein Description:** Human GHRHR-Strep full length protein-synthetic nanodisc

**Formulation:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.

**Protein Pathways:** GPCRDB Class B Secretin-like, GPCRDB Other, Cancer.

**Protein Families:** GPCR, Transmembrane, Druggable Genome.

**Usage:** Research use only

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.