

RPCB0025

Product Information

Product SKU:	RPCB0025	Gene ID:	269	Size:	10µg
Tag:	C-hFC	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q16671
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: AMHR2, On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human AMHR2/MISR2 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.01 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.