

Biotinylated Recombinant Human C1qR/CD93 Protein

RPCB0070



Product Information

Product SKU:	RPCB0070	Gene ID:	22918	Size:	100µg
Tag:	C-His&Avi	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q9NPY3
Purity:	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

Protein Information

Background: CD93 has been shown critical roles in inflammatory and immune diseases. CD93 silencing increased IL-6 and TSLP, but not IL-33 levels in culture supernatants. HDM-induced asthma mice showed significant airway hyperresponsiveness and inflammation with Th2 cytokine activation, along with decreased CD93 expression in bronchial epithelial cells and lung homogenates but increased serum CD93 levels.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human C1qR/CD93 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.