

Biotinylated Recombinant Human Semaphorin-4D/SEMA4D/CD100 Protein

RPCB0073



Product Information

Product SKU:	RPCB0073	Gene ID:	10507	Size:	100µg
Tag:	C-His&Avi	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q92854-1
Purity:	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

Protein Information

Background: Semaphorin 4D (Sema4D) is a multifunctional protein widely expressed in an organism that plays an important role in the control of many physiological and pathological processes, including immunoregulation, neurogenesis, angiogenesis, and tumor progression.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human Semaphorin-4D/SEMA4D/CD100 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.