

Biotinylated Recombinant Human TGF-beta receptor type-1/TGFR-1/ALK-5 Protein

RPCB0083



Product Information

Product SKU: RPCB0083 **Gene ID:** 7046 **Size:** 100µg
Tag: C-His&Avi **Reactivity:** Human

Additional Information

Expression Host: - **Swissprot:** P36897-1
Purity: > 95% by Tris-Bis PAGE

Protein Information

Background: transforming growth factor beta receptor 1 (TGFB1), a key stimulator of tumor proliferation and metastasis, was a direct target of miR-98-5p. miR-98-5p overexpression resulted in the downregulation of TGFB1 and the suppression of the viability, proliferation, migration and invasion of A549 and H1299 cells.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human TGF-beta receptor type-1/TGFR-1/ALK-5 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.