

# Biotinylated Recombinant Mouse TNFRSF17/BCMA/CD269 Protein

RPCB0089



## Product Information

**Product SKU:** RPCB0089      **Gene ID:** 21935      **Size:** 100µg  
**Tag:** C-His&Avi      **Reactivity:** Mouse

## Additional Information

**Expression Host:** -      **Swissprot:** O88472-1  
**Purity:** > 95% by Tris-Bis PAGE and HPLC

## Protein Information

**Background:** B-cell maturation antigen (BCMA or BCM), also known as tumor necrosis factor receptor superfamily member 17 (TNFRSF17), is a protein that in humans is encoded by the TNFRSF17 gene. TNFRSF17 is a cell surface receptor of the TNF receptor superfamily which recognizes B-cell activating factor (BAFF).

**Protein Description:** High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Mouse TNFRSF17/BCMA/CD269 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Normally 8% trehalose is added as protectant before lyophilization. Contact us for customized product form or formulation.

**Storage:** Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.