

# Biotinylated Recombinant Human Angiopoietin-like 2/ANGPTL2 (260-493) Protein

RPCB0095



## Product Information

<b>Product SKU:</b>	RPCB0095	<b>Gene ID:</b>	23452	<b>Size:</b>	100µg
<b>Tag:</b>	C-His&Avi	<b>Reactivity:</b>	Human		

## Additional Information

<b>Expression Host:</b>	-	<b>Swissprot:</b>	Q9UKU9-1
<b>Purity:</b>	> 95% by Tris-Bis PAGE and HPLC		

## Protein Information

<b>Background:</b>	Angiopoietin-like 2 (ANGPTL2) is a proinflammatory protein belonging to the angiopoietin-like family. ANGPTL2 is secreted and detected in the systemic circulation. Different observational clinical studies reported that circulating levels of ANGPTL2 increase significantly in various chronic inflammatory diseases and showed associations between ANGPTL2 levels and diagnosis and/or prognosis of cardiovascular diseases, diabetes, chronic kidney disease, and various types of cancers.
<b>Protein Description:</b>	High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human Angiopoietin-like 2/ANGPTL2 (260-493) Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.
<b>Endotoxin:</b>	< 1 EU/µg of the protein by LAL method.
<b>Formulation:</b>	Lyophilized from a 0.22 µm filtered solution of PBS, 200 mM L-arginine, pH 7.4. Normally 8% trehalose is added as protectant before lyophilization. Contact us for customized product form or formulation.
<b>Storage:</b>	Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.