

Product Information

Product SKU:	RPCB0129	Gene ID:	81607	Size:	100µg
Tag:	C-His&Avi	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q96NY8
Purity:	> 95% by Tris-Bis PAGE and HPLC		

Protein Information

Background: Nectin-4 seems to be involved in cell adhesion through trans-homophilic and - heterophilic interactions, the latter including specifically interactions with NECTIN1. Does not act as receptor for alpha-herpesvirus entry into cells.(Microbial infection) Acts as a receptor for measles virus.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human Nectin-4/PVRL4 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Normally 8% trehalose is added as protectant before lyophilization. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.