

RPCB0201

Product Information

Product SKU:	RPCB0201	Gene ID:	5176	Size:	20µg
Tag:	No tag	Reactivity:	Human		

Additional Information

Expression Host:	Baculovirus-Infected Sf9 Cells	Swissprot:	P36955
Purity:	> 97% by SDS-PAGE.		

Protein Information

Background: The protein is a member of the serpin family, although it does not display the serine protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in neuronal differentiation in retinoblastoma cells.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Serpin F1/PEDF Protein, tested reactivity in Baculovirus-Infected Sf9 Cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.