

Recombinant Human Alpha-2-MRAP/LRPAP1 Protein

RPCB0214



Product Information

Product SKU:	RPCB0214	Gene ID:	4043	Size:	10µg
Tag:	C-His	Reactivity:	Human		

Additional Information

Expression Host:	E. coli	Swissprot:	P30533
Purity:	> 87% by SDS-PAGE.		

Protein Information

Background: This protein that interacts with the low density lipoprotein (LDL) receptor-related protein and facilitates its proper folding and localization by preventing the binding of ligands.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Alpha-2-MRAP/LRPAP1 Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of 20mM Tris, 150mM NaCl, 5% glycerol, pH 8.0. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.