

RPCB0254

Product Information

Product SKU: RPCB0254	Gene ID: 7037	Size: 10µg
Tag: C-His	Reactivity: Human	

Additional Information

Expression Host: HEK293 cells	Swissprot: P02786
Purity: > 97% by SDS-PAGE.	

Protein Information

Background: This protein is a cell surface receptor necessary for cellular iron uptake by the process of receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Multiple alternatively spliced variants have been identified.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human sTfR/CD71 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.