

RPCB0428

Product Information

Product SKU:	RPCB0428	Gene ID:	3491	Size:	10µg
Tag:	C-hFc	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	O00622
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: The secreted protein encoded by this gene is growth factor-inducible and promotes the adhesion of endothelial cells. The encoded protein interacts with several integrins and with heparan sulfate proteoglycan. This protein also plays a role in cell proliferation, differentiation, angiogenesis, apoptosis, and extracellular matrix formation.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Cyr61/CCN1 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.1 EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.