

RPCB0440

Product Information

Product SKU:	RPCB0440	Gene ID:	54567	Size:	10µg
Tag:	C-His	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	Q9NR61
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: This protein is a homolog of the Drosophila delta family. The delta family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human DLL4 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.