

RPCB0465

---

## Product Information

<b>Product SKU:</b>	RPCB0465	<b>Gene ID:</b>	4283	<b>Size:</b>	10µg
<b>Tag:</b>	No tag	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Expression Host:</b>	E. coli	<b>Swissprot:</b>	Q07325
<b>Purity:</b>	> 90% by SDS-PAGE.		

---

## Protein Information

**Background:** This antimicrobial protein belongs a protein thought to be involved in T cell trafficking. The encoded protein binds to C-X-C motif chemokine 3 and is a chemoattractant for lymphocytes but not for neutrophils.

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human CXCL9/MIG Protein , tested reactivity inE. coliand has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 0.1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.