

RPCB0491

Product Information

Product SKU: RPCB0491	Gene ID: 5817	Size: 100µg
Tag: C-mFc	Reactivity: Human	

Additional Information

Expression Host: HEK293 cells	Swissprot: P15151-1
Purity: > 95% by SDS-PAGE ; > 95% by HPLC	

Protein Information

Background: CD155 is a cell surface adhesion molecule functioning in tumor cell migration, invasion, and metastasis, and not surprisingly, is also designated as a common tumor-associated antigen. CD155 is also recognized by NK cells to induce their cytotoxicity. CD155 is also commonly referred to as the "poliovirus receptor," or PVR.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human PVR/CD155 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.