

RPCB0517

Product Information

Product SKU:	RPCB0517	Gene ID:	4680	Size:	100µg
Tag:	C-His&Avi	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P40199
Purity:	> 95% by SDS-PAGE ; > 95% by HPLC		

Protein Information

Background: Carcinoembryonic antigen-related cell adhesion molecule 6 (CEACAM6) belongs to the human carcino-embryonic antigen (CEA) family. Numerous lines of studies have indicated that altered expression of CEACAM6 may have a role in carcinogenesis and development.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human CEACAM6/CD66c Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.